

Phenyl

Mechanisms of Separation

π - π interactions

Strength of Interaction

Strong

Dipole-dipole interactions

Moderate

Hydrophobic binding interactions

Moderate

Target Analytes

Analytes with π bonding

Analytes with electron delocalization and electron-withdrawing groups, such as halogens, nitro groups, ketones, esters, and acids

Analytes with proton donor groups

Analytes with different dipole moments

Recommended Applications

Stereoisomers

Steroids

Taxanes

Substituted aromatics


Highly aqueous conditions

FIGURE 27: Comparison of Selectivity differences between C18 and Phenyl phases


Column: 4.6 x 250 mm, 5.0 μ m
Mobile Phase: 20% 0.025M KH_2PO_4 , pH 6.0
80% MeOH
Flow Rate: 1.0 mL/min
Temperature: 22 $^\circ\text{C}$

1. Norephedrine
2. Nortriptyline
3. Toluene
4. Imipramine
5. Amitriptyline

ACE C18



ACE Phenyl



These comparison chromatograms show the substantial differences in selectivity between a C18 and a phenyl phase. Peaks 2 and 3 have reversed elution order on the phenyl phase compared to the C18 phase