

## Precolumn Filters

Placed immediately before the column, precolumn filters trap sample particulates. However, for samples that may irreversibly adsorb onto the column, guard cartridges are preferred.

- Protect columns from particulates
- Compatible with all HPLC columns
- Ultra low dead volume

### ColumnSaver precolumn filter

The ColumnSaver fingertight precolumn filter is universally compatible with all column manufacturers' endfittings. Both  $2\mu\text{m}$  and  $0.5\mu\text{m}$  versions are available, for the protection of  $5\mu\text{m}$  and  $3\mu\text{m}$  columns respectively. The ColumnSaver precolumn filter is leakproof to over 5000psi.



Part No.	Description	Qty
HI-685	$2\mu\text{m}$ ColumnSaver	10/pk
HI-686	$0.5\mu\text{m}$ ColumnSaver	10/pk

### Stand-alone precolumn filter

Model HI-704 comprises a 'stand-alone' holder with a  $2\mu\text{m}$  PEEK encapsulated stainless steel frit.



### Precolumn filters

Part No.	Description	Qty
HI-704	$2\mu\text{m}$ Precolumn filter assembly (complete with frit)	1/pk
HI-101	Replacement $2\mu\text{m}$ frits	10/pk
HI-102	Replacement $0.5\mu\text{m}$ frits	10/pk

## HPLC Accessories