


## Conditions

Column: ACE 3 C18-AR  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-119-1546  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: 0.1% formic acid in MeOH

Time (mins)	%B
0	52
28	74
33	74
38	52
48	52

Flow Rate: 1 mL/min  
Injection: 5 µL  
Temperature: 40 °C  
Detection: UV, 254 nm



1. Bendroflumethiazide
2. Ketoprofen
3. Naproxen
4. Sulindac
5. Ibuprofen
6. Diclofenac
7. Indomethacin
8. Meclofenamic acid
9. Mefenamic acid